







DATE	SYM	REVISION RECORD	ECO	DWN	ENG	MFG
2025-02-12	01	INITIAL RELEASE		EB	JS	FD



ITEM NO.	PART NUMBER	DESCRIPTION	QTY.
1	FFH-0202-VB0	Internal Ferrule 1/16 SS	4
2	FNN-S202-VA8	Internal Nut 1/16" - Short-	4
3	ROT-00015-01-A	MRV ISL 4-PORT 0.1µL ROTOR ASSEMBLY	1
4	ROT-00018-01-A	ISL 4-Port MRV Stator With 1/16 Port and 0.015" Bore	1
5	RV00-0198-1000	RV/MRV Short Preload Assembly	1

UNLESS OTHERWISE SPECIFIED: - ALL DIMENSIONS ARE IN INCHES - SURFACE FINISH IS IN Ra MICROINCHES (µIN) - UNSPECIFIED SURFACES MUST BE 64 Ra µIN OR BETTER - REMOVE ALL BURRS & SHARP EDGES .003-.010 - TOLERANCES ARE: - XX ±.010 FRACTIONAL ±1/64" - XXX ±.003 ANGULAR ±5° - XXXX ±.001 - UNSPECIFIED INSIDE CORNER RADIUS 0.010 MAX - REFER TO 3D MODEL FOR DIMENSIONS NOT SHOWN, TOLERANCE XX APPLIES - DO NOT SCALE DRAWING	MATERIAL:		 2659, boul. du Parc Technologique, Québec Québec, CANADA G1P 4S5 Tel: 418-338-0004 Fax: 418-338-2500 Email: info@afproducts.ca Web: www.afproducts.ca			
	FINISH:					
	X		DESCRIPTION:			
	ACCORDANCE WITH LATEST ASME Y14.5		MRV-4 ISL 0.1µL 1000psi			
PROPRIETARY INFORMATION THIS DRAWING IS THE PROPERTY OF APN-AFP INC. IT IS RETURNABLE UPON REQUEST AND MUST NOT BE REPRODUCED OR THE INFORMATION CONTAINED HEREIN TRANSMITTED TO ANY UNAUTHORIZED PERSON OR ORGANIZATION WITHOUT THE WRITTEN CONSENT OF APN-AFP INC.			UNIT	SIZE	DWG. NO.	REV.
			IN	C	MRV04-16R2-0023	01
			SCALE	3:2	SHEET	1 OF 2