







Date	Ver-Rev.	Description - ECR	Design by
2025-07-22	01	INITIAL RELEASE	Etienne Boisclair



SPECIFICATIONS:

PRODUCT P/N	MRV04-16R2-0026-CP
VALVE P/N	MRV04-16R2-0026
VALVE TYPE	HIGH PERFORMANCE GRADE, MEDIUM TEMPERATURE
NUMBER OF PORTS	4
PORT SIZE	0.025" [0.64 mm]
PORT CONNECTION	1/16" [1.6 mm]
INTERNAL SAMPLING LOOP SIZE	2µL
MAX. WORKING PRESSURE	1000 psig [6895 kPa]
MAX. INITIAL LEAK RATE	1.6E-10 Atm-cc/s He
MAX. INITIAL PRESSURE DROP	2 kPa @ 75 sccm N ₂ , 1 kPa between any section
STAND-OFF LENGTH	CLOSE MOUNT
ACTUATOR PRESSURE INLET	40 psig [275 kPa] min., 80 psig [550 kPa] max.
VALVE MAX. WORKING TEMPERATURE	437°F [225°C]
ACTUATOR MAX. WORKING TEMP.	257°F [125°C]
MATERIAL	NITRONIC 60 VALVE STATOR R2 ROTOR TYPE (MODIFIED PEEK)

 2659, boul. du Parc Technologique Québec, Québec, CANADA, G1P 4S5 418-338-0004	Part MRV04-16R2-0026-CP		Material	Finish	
	Description Closemount MPA 2µL 1000psi				
	Rev. 01	Status Released	Design by Etienne Boisclair		Approved by à valider avec Josef
	Unit INCH	Format C	Scale 1:1		Sheet 1 of 2
	ACCORDANCE WITH LATEST ASME Y14.5				
	PROPRIETARY INFORMATION: THIS DRAWING IS THE PROPERTY OF APN-AFP INC. IT IS RETURNABLE UPON REQUEST AND MUST NOT BE REPRODUCED OR THE INFORMATION CONTAINED HEREIN TRANSMITTED TO ANY UNAUTHORIZED PERSON OR ORGANIZATION WITHOUT THE WRITTEN CONSENT OF APN-AFP INC.				